**Intended use**

Easicult M slides are intended for monitoring of fungal contamination in various industrial environments. The test can be performed on-site, or the slides can be used as convenient transport media for samples.

The slide is covered on both sides with a growth medium which supports growth of both moulds and yeasts. The main significance of the test is that elevation of fungal density can be detected. As no universally applicable limits for fungal density exist, limit values have to be determined by experience.

**Contents of the kit**

Easicult M	Cat. No. 67686
Test slides	10 pcs
Labels	10 pcs
Instructions for use	1 pc

**Typical formulation of growth medium**

Malt agar	Chloramphenicol
Yeast extract	Gentamycin sulphate
Dextrose	Agar-agar
Lactic acid	Water

**Warnings and precautions**

Do not use the product beyond the expiry date marked on the kit.

Do not touch the unused growth media.

Do not use the slide if you notice

- discoloration or dehydration of the growth medium
- detachment of the growth media from the plastic slide
- evidence of microbial growth.

Because any growth on the Easicult M slide may be pathogenic, do not touch the used slide.

**Storage**

Store Easicult M at room temperature (approx. 20°C/68°F) protected from draught, temperature fluctuations and light sources. Avoid storage near heat-generating appliances. Do not allow to freeze. The expiry date (year-month-date) is marked on the box and on the cap of each slide.

**Sampling and procedure (Fig. 1–5)**

To avoid contamination, the growth medium should not come into contact with any other material than the one to be tested. On the other hand, it is important that the growth medium makes full contact with the material to be tested.

**Viscous fluids and fluids with high microbial content**

If the viscosity or microbial content of the sample is high, the sample should be diluted. For dilution, put 100 or 1000 ml of drinkable tap water in a clean, well-rinsed and dried bottle with a cap. The bacterial content of the water for dilution should not exceed 100 CFU/ml. Before filling the bottle, let the water run for 5 minutes or boil it for 15 minutes and then let it cool. Using a clean (disposable) pipette, add 1 ml of sample to the bottle. Close and mix thoroughly by shaking the bottle about 30 times. Dip the slide in the dilution and proceed as described for fluid samples.

**Fluid samples**

- 1 Unscrew the tube and withdraw the slide without touching the agar surfaces.
- 2 Dip the slide in the liquid. Alternatively, wet the slide under a running stream of the liquid or spray the slide. If the liquid is under pressure, the slide must be handled carefully to avoid unfastening of the agar. If the sample is in a container, mix the contents and dip the slide in it.  
Both agar sides should get completely wet. The slide must be in contact with the fluid for 5 to 10 seconds.
- 3 Allow excess fluid to drain off the slide. Blot the last drops from the lower end of the slide on absorbent paper.
- 4 After sampling, screw the slide tightly back into the tube. Fill in the label and affix it to the tube.
- 5 Incubate the slide at 27...30°C (80...86°F) for two days. Some slow-growing organisms may not be visible after two days' incubation. If there is no growth within two days, further two days of incubation are recommended.

**Interpretation of results (Fig. 6)**

Cautiously remove the slide from its tube after incubation and determine the density (number of colonies) and characteristics of growth by comparing the growth with the model chart.

If the sample was diluted, the dilution factor must be taken into account in the evaluation. For example, if a dilution of 1 + 100 ml (1 ml of sample in 100 ml of water) shows for yeasts a density of 10<sup>6</sup> CFU/ml, the actual result for the sample is 10<sup>8</sup> CFU/ml (CFU = number of colony forming units).

Growth appearing on the slide may consist exclusively of either yeasts or moulds, or mixed growth of both. Mould colonies are soft and fluffy and usually whitish, green or black in colour. Yeasts usually grow in dome-shaped colonies but may sometimes be flat and dry. Yeast colonies are often red or white. Since mould colonies may originate from fragments of mycelium or from individual spores, the model chart is not quantitative. The chart indicates whether the contamination is slight (+), moderate (++) or heavy (+++). Colonies can be removed from the slides and examined under microscope if needed. Fungal contamination is sometimes apparent to the naked eye as a coating on the fluid surface.

**Limitations of the method**

The reliable lower detection limit for yeasts is 10<sup>3</sup> CFU/ml.

**Disposal**

- Dispose of contents according to national and local law.
- All patient samples and used components should be handled and disposed of as potentially infectious material.
- Materials of the components:  
Paper: Instructions for use, patient labels  
Cardboard: Kit box  
Plastic: Tubes, caps and dipslides
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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Easicult® is a registered trademark of Orion Diagnostica Oy.